



Proud exclusive distributors of Biobase equipment & accessories in Southern Africa

Tel: (011) 763-8820/8405    FAX: (011) 763-2769    e-mail: [sales@cleanroom.co.za](mailto:sales@cleanroom.co.za)/[crm@corpdial.co.za](mailto:crm@corpdial.co.za)  
Emergency number: 0829223565 (Cornus)    website: [www.cleanroom.co.za](http://www.cleanroom.co.za)

## Air Curtain Tester



Introductory Price

R950.00

Normal Price R1150.00

### **Introduction:**

This equipment is used to test the airflow patterns within a Microbiological safety cabinet and/or a laminar flow unit. This confirmation of airflow patterns ensures the unit's optimal working can be checked in-house and thus ensures product & operator safety.

It works on the same principle as the smoke test tubes but only needs water & electricity to work. It is more durable, user friendly & cost effective.

### **Air Curtain Tester**

Power supply: 220 AC 50HZ

Input power: 30w Fuse: T2A L (  $\phi$  5×20)

Work style : Intermittent working≤15 min

Maximum atomization velocity: 1.5 ml/min

Maximum volume: 35ml

### **Working environment:**

Environment temperature: 10-40°C

Relative humidity: ≤70% RH

Air pressure range: 860-1060 kPa

### **Storage environment:**

Environment temperature: 10-40°C

Relative humidity: ≤95% RH

Air pressure range: 500-1060 kPa

**For sales enquiries contact Elliot 082 446 2258 or Cornus 082 922 3565**



Proud exclusive distributors of Biobase equipment & accessories in Southern Africa

Tel: (011) 763-8820/8405 FAX: (011) 763-2769

e-mail: [sales@cleanroom.co.za](mailto:sales@cleanroom.co.za)/[crm@corpdial.co.za](mailto:crm@corpdial.co.za)

Emergency number: 0829223565 (Cornus)

website: [www.cleanroom.co.za](http://www.cleanroom.co.za)

## Formalin Fumigation Sterilizer and Ammonium Hydrogen Carbonate Neutralizer



### Introduction

These two equipments are used for completely decontamination of the microbiological safety cabinet before changing filters and/or major maintenance..

#### Formalin Fumigation Sterilizer:

Model: FX-100

Voltage/Frequency: 220V/50 Hz

Power: 180w

Volume: 100 ml

It used for sterilization of the safety cabinet before changing filter

#### Ammonium Hydrogen Carbonate Neutralizer:

Model: TZ-100

Voltage/Frequency: 220V/50 Hz

Power: 180 w

Volume: 100 ml

It used for toxicity elimination of the Formalin

**For sales enquiries contact Elliot 082 446 2258 or Cornus 082 922 3565**